



Labosport US Inc • 1806 S. Dixie Highway, Suite G • 30720 • Dalton • Georgia

## SAMPLING INSTRUCTIONS

Please follow the instructions below **to collect a representative infill sample from your field:**


- Use the diagrams below to collect infill from different locations of your field.
- Collect a **total weight of approx. 20 oz of infill** (for instance, take 3-4 oz from each locations)
  - Use an exam glove to avoid any external contamination of the infill (food residues for instance).
- Ensure that the rubber infill is free from sand, as much as possible.
- Place your infill sample in a closed transparent clean and single use plastic bag (ex: ZIPLOC Bag).
- Put the bag in a box and ship it using the tag provided by our customer service.



### Football



### Baseball



### Soccer & others



Within 48h you will received a shipment tag to ship your sample.

For any question please contact Labosport USA customer service at [customerserviceusa@labosport.com](mailto:customerserviceusa@labosport.com) or call at 706-529-9474